

# iLLUNA PLLA

200mg x lvial 500mg x lvial



iLLUMA PLLA rejuvenates the skin by stimulating collagen production, resulting in natural volume & firm, healthy skin

# Composed of porous micro-particles

200mg x vial 500mg x vial



Promotes tissue regeneration and excellent collagen production with minimal application

Maintains a natural volumizing effect for a long time

Consistent particle shape enables smooth injection without needle clogging

It stimulates collagen production within the face and body, thereby preserving the innate suppleness and elasticity of the skin



#### **PRODUCT FEATURES**

**VOLUME: 200mg / 500mg DILUTION VOLUME: 8ml / 20ml CANNULA: ILLUMA PLLA 200 Use** cannula 27-28g iLLUMA PLLA 500 Use cannula 23-**25g** 

> MAIN INGREDIENTS PLLA (Poly-L-Lactic Acid) **CMC Compound Amount**

> NUMBER OF TREATMENTS 1-3

## Advantages of iLLUMA PLLA

The focus on natural and safe results has led to increased consumer interest, attracting and retaining patients.

### Safety

The product, utilizing the NST method, is eco - friendly and solventfree, ensuring safe and stable facial procedures.

#### **Natural Outcomes**

Made from natural corn-derived polymer, the product provides a voluminous effect and stimulates collagen regeneration.

# **Time-Efficiency**

The significantly reduced dilution time cuts preparation to just 20 minutes. enhancing satisfaction and efficiency for both users and patients

